

## **Storage of data on the Bioimg server**

### **a. Exporting the data:**

1. Create in "D:\UserData\" a folder with your exact User name (in small letters!)
2. Export your data from BD FACSDiva in FCS file format or in Experiment format (depending on requirements) into the folder created.
3. Wait for the data to be backed up by the Bioimg server (up to two hours)

### **B.1. Accessing the data: (PC)**

1. Right click on "My Computer"
2. Choose "Map Network Drive"
3. Choose "Drive: Z:"
4. Folder: "\\bioimg\User name"
5. Click on "Finish"
6. Wait a few seconds until "**Connect to bioimg**" window appears
7. Write your **User Name**
8. Write the Password **User Name1**
9. The folder of your user name on bioimg data will open.  
Copy the folder required by dragging it.

### **B.2. Accessing the data: (MAC)**

1. Go to "Go>Connect to server"
2. Write: "smb://bioimg"
3. Wait a few seconds until **Connect to bioimg** window appears
4. Write your **User Name**
5. Write the Password **User Name1**

### **C. Disconnect from bioimg**

1. Open "My Computer"
2. Right click on your User Name on 'BioImaging Data (bioimg)' (Z:) and disconnect
3. On the MAC put the connected server in the trash

### **D. Accessing and managing the data on the BioImaging Storage Server**

1. Open the chrome internet browser and go to "bioimg.weizmann.ac.il"
2. Logon with your User name and your regular Weizmann password
3. Use "annotate Raw Data" to add experiment to the database by assigning it to a project.
4. Use "Frequently asked Questions" for help